

## Easy to Use Forehead Thermometer

### Safe, Easy, Quick Temperature Monitoring

The **NEW** Wahl Fever Alert Forehead Thermometer is a liquid crystal forehead thermometer which provides an accurate and easy to read dual scale temperature reading by changing color to indicate an elevated skin temperature. You can obtain an accurate temperature in approximately 15 seconds.

Simply watch the liquid crystal display until the colors stop changing. Since you are able to take a subject's temperatures in a non-invasive way, the forehead thermometers are ideal for use with anyone, allowing for non-contact fever screening.

This disposable thermometer is an FDA Class II registered medical device with a reusable and repositionable medical grade adhesive. They are a cost and time-effective solution for a quick check of an individual's temperature. These thermometers are reusable (up to 5 times by the same person) allowing you to monitor a subject's temperature throughout the day or week and still maintain an appropriate social distance.

As a medical device, the forehead label can also be used as a tool for verification if a thermal imaging or infrared system indicate an elevated skin temperature in an individual. A physician should always be consulted if a fever is detected.

#### Directions for Use:

Forehead thermometers should be used in-doors at room temperature (approximately 68° to 86°F). Avoid direct sunlight or strong heat.

Do not take the temperature of an individual if they have been eating, drinking or exercising within the last 30 minutes due to the possibility of temperature fluctuation.

1. Remove backing.
2. Hold thermometer at both ends, pressing firmly on dry forehead for approximately 15 seconds.
3. Read Thermometer while on forehead. Green indicates the correct temperature.

If you do not see Green, Blue or Tan will appear:

Blue: add 1°F/0.5°C / Tan: subtract 1°F/0.5°C.

4. Reusable up to 5 times.

Always consult your physician if a fever is detected.

**To be administered by adults only. Keep out of reach of children.**

**See how  
it works!**



### Specifications

Temperature points	95°F/35°C 96.8°F/36°C 98.6.4°F/ 37°C 100.4°F/38°C 102.2°F/39°C 104°F/40°C
Accuracy	± 1.8°F / 1°C
Classification	FDA Class II Medical Device Class I in accordance with EU 2017/745. ISO 9001:2015
Label Type	Reversible, Liquid Crystal based. Reusable with repositionable medical grade adhesive.
Dimensions	2.825 in x .44 in (7.1755 cm x 1.12 cm)
CE Approval	Yes
Shelf Life	2 years



Available in packages of 10/20/100 with reclosable zipper top. (Package includes hang hole which fits standard peg hooks for retail).

Also available in rolls of 5,000.

Labels are perforated for easy separation.

**Important Note:** The Forehead Thermometer is a medical device that measures skin temperature. It cannot be used to detect or diagnose COVID-19/Corona virus.

*Specifications subject to change without notice.*

June  
2020

Palmer Wahl Instruments, Inc. • Toll Free: 800-421-2853 • Phone: 828-658-3131 • www.palmerwahl.com

